


EN	OVATION HYBRID NIGHT SPLINT®
DE	OVATION HYBRID-NACHTSCHIENE®
ES	FÉRULA NOCTURNA HÍBRIDA OVATION HYBRID NIGHT SPLINT®
FR	ATTELLE DE NUIT HYBRIDE OVATION®
IT	STECCA NOTTURNA IBRIDA®
PL	ORTEZA NA OPADAJĄCĄ STOPEŃ OVATION HYBRID NIGHT SPLINT®

ENGLISH

BEFORE USING THE DEVICE, PLEASE READ THE FOLLOWING INSTRUCTIONS COMPLETELY AND CAREFULLY. CORRECT APPLICATION IS VITAL TO THE PROPER FUNCTIONING OF THE DEVICE.

INTENDED USE/INDICATIONS: The Hybrid Night Splint® is intended for use as a support for the ankle and toe region as well as fixation of foot ankle and stretching of the plantar fascia during inactivity. The product is indicated for plantar fasciitis, Achilles tendonitis, drop foot and post-static pain.

The Hybrid Night Splint® is intended for short term single patient use and not intended to undergo regular service cycles.

This product was designed to withstand normal use conditions for orthopedic support, post injury rehab or prevention of further injury. Like all products that are worn by individuals, the extended life of such products is determined by the level of activity and type of activity of the individual.

CONTRAINDICATIONS: None known. In every case, the final decision whether to use this Ovation Medical® device rests with the attending physician/therapist. He/she should assess prior to use of the device, if it is suitable for a given application, taking into account the condition of the patient.

WARNINGS AND PRECAUTIONS: In the event of discomfort or numbness of any kind when using this product, consult your health care provider immediately. While non-slip allows for safe ambulation, the Hybrid Night Splint® is not intended for continuous day-in/day-out walking. Please consult your health care provider for specific instructions. This product is to be used under the supervision of a medical professional.

REPORTING PROBLEMS: The user should report any suspected serious incident related to the use of this Ovation Medical® device by informing the manufacturer and the competent authority of the state or the country in which the incident occurred.

APPLICATION INSTRUCTIONS:

- To begin, unfasten all straps.
- Use either method (4a or 4b) to obtain desired tension/position.
- Use position indicator to measure progress.

CARE AND MAINTENANCE: Hand wash in cold water using mild soap, do not bleach, line dry only, avoid direct heat or sun. Store in a cool dry place.

DEVICE DISPOSAL: Ovation Medical® products should be properly disposed of by a qualified recycler when the post-operative use period is completed or when the device is not longer required by the patient.

SYMBOLOLOGY:

- Consult instructions for use and accompanying documents.
- Manufactured by. Note: can be used in combination with date of manufacture.
- Date of manufacture accompanied by a date.
- Single patient multiple use.
- Device sold in non-sterile condition.
- Not made with Natural Rubber Latex
- Peel here.
- Lot number.
- Reference number or product number.
- Quantity contained in package.
- UDI Unique Device Identifier Number.
- Product intended to be used as a medical device.
- Authorized representative (REP) in the European Community (EC).
- Conform to the European Medical Device Regulation (EU) 2017/745 for class I medical devices

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