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Dear customer,  
Thank you for choosing a Bauerfeind product.

We work to improve the medical effectiveness of our products every day – because your health is very important to us. Please read and observe these instructions for use carefully. If you have any questions, please contact your doctor or medical retailer.

### Intended purpose

SecuTec Genu is a medical product. It is a stabilizing orthosis that limits the range of motion in the case of complex knee injuries.

### Indications

- Rupture of the anterior and / or posterior cruciate ligament (ACL / PCL)
- After ligament surgery / ligament reconstruction
- Severe and / or complex instabilities (traumatic, degenerative, e.g. "unhappy triad")
- Collateral ligament injuries
- To restrict the knee joint's range of motion (e.g. after meniscus resection or meniscal implantation)
- Patellar fracture (non-surgical and post-operative)

### Risks of using this product

- Once you have been prescribed the SecuTec Genu, only use it according to your indication and any additional instructions given by a medical specialist<sup>2</sup>. Please consult a specialist or your physician if you want to use it together with other products. Do not modify the product yourself in any way; otherwise it may not help as expected or it may be harmful to your health. In these cases, any warranty or liability claims will be excluded.
- Avoid contact with ointments, lotions, or substances containing grease or acids as well as extreme temperatures.
- Any aids applied externally to the body may, if tightened excessively, lead to local pressure points or, in rare cases, constrict the underlying blood vessels or nerves.
- If you notice any changes or an increase in symptoms while wearing the product, stop any further use and contact your doctor.
- Should any problems with the fit of the orthosis arise, for instance if the shape of your leg changes (e.g. due to an increase or reduction in muscle), or should the product have a functional defect, please immediately contact your medical supply retailer.
- No side effects that affect the entire body have been reported to date. This assumes the orthosis is used / fitted correctly.

### Contraindications

No clinically significant adverse reactions have been reported to date.

If you have any of the following conditions, the product should only be used after consultation with your physician:

- Skin disorders / injuries to the relevant part of the body, particularly if inflammation is present. Likewise, any raised scars with swelling, redness, and excessive heat build-up
- Varicose veins (varicosis)
- Impaired sensation and circulatory disorders of the legs / feet (e.g. diabetes mellitus).
- Impaired lymph drainage – including soft tissue swellings of uncertain origin located away from the fitted aid.

### Application instructions

#### Size determination

The size is determined by trained specialists according to the requirements. You can find a size chart on the product packaging or on our homepage [www.bauerfeind.com](http://www.bauerfeind.com).

#### Procedure for the first fitting of the SecuTec Genu orthosis

#### Adjustment options for the hinge (to be performed by a trained specialist only)

- The hinge cover can be opened via a screw fastening. This also serves to fasten the movement-limiting wedges (extension – front / flexion – rear) and protects the hinge mechanism against dirt.
- Opening the screw fastening (quarter turn until the arrow points to ) allows the hinge protection to be removed and the desired limiting wedges to be changed.
- The desired extension and flexion wedges are selected from the enclosed frame and are used ①②. The following limitations are possible: Extension: 0°, 10°, 20°, 30°, 45° Flexion: 0°, 10°, 20°, 30°, 45°, 60°, 75°, 90°
- Adjustments must be carried out in the same way on the inside and the outside.
- Then, click the hinge cover onto the hinge, check whether it is correctly positioned and firmly attached to the hinge, and close the screw fastening  ③. Make the adjustment to both hinges.

#### Static fitting (only by a trained specialist)

- Perform the initial adjustment with the straps open.
  - Bend the knee to be treated to a 30° angle and place the SecuTec Genu on it so that the hinge pivot point of the orthosis hinge is aligned with the compromise axis of rotation as defined by Nietert ④.
  - Check the fit of the frame parts when they are laid flat to make sure their circumference and longitudinal direction are correct and adjust them if necessary.
  - A bending iron with protective caps should be used for any extensive bending work.
  - Select the required thickness of knee joint cushions and insert them (there are two thicknesses to choose from).
- Frontal version:

- If necessary, the shin cushion provided can be attached in order to redistribute pressure in the event of pressure pain in the shin region or swelling around the tibial tubercle (Tuberositas tibiae) ⑤.
- With the "frontal" version, you decide on the calf strap height that provides the optimum fit for the patient. To do so, remove the calf strap from the frame and move the stopper in the slot. For instance, if the stopper is pushed upwards, pull the strap back through the slot and hook the quick-fit fastener into the lower loop. The strap is now adjusted for a low-lying calf. Now hook in the remaining quick-fit fasteners.

#### Putting on the SecuTec Genu (for the patient)

- Put the orthosis on against your skin. Wearing the support over pants can cause slipping and impair functionality.
- Put on the "dorsal" version of the orthosis by holding it in front of you, placing your foot down through the frame, and then pulling the frame up to your knee. Put on the "frontal" version of the orthosis by placing it on your extended leg.
- The leg is bent at a gentle angle (approximately 30°). The position of the kneecap is then determined through palpation.
- The center of the hinge (pivot point of the orthosis hinges) is positioned at approximately the same height as the center of the kneecap ⑥.
- Start by pulling the strap tight directly beneath the knee, then proceed to the straps closer to the joint above the knee. – When closing the straps, check whether the orthosis has rotated laterally around the leg and correct this if necessary.
- Position the front strap with the protective cushion for the tibial tubercle in such a way that the strap can be closed without exerting any pressure. After all the straps have been closed, the orthosis must exert an even pressure when the muscles are tensed.
- Once the orthosis has been put on, check once more that it is in the correct position.

#### Taking off the SecuTec Genu

- Grip the handle of the quick-fit fastener and unhook it from the thigh frame.
- Then lift the upper orthosis frame off the thigh. During this step, the knee joint cushions stay in contact with the knee joint.
- Now remove the straps from the lower leg frame of the orthosis.
- The straps without quick-fit fasteners can remain closed.
- After all the straps have been removed, the orthosis can be taken off the knee joint from the front.
- After all the straps have been removed from the "frontal" version, the orthosis can be taken off the knee joint from the front. For the "dorsal" version, push the orthosis down toward your foot and step out of the frame.

#### Cleaning instructions

#### Never expose the orthosis to direct heat (e.g. central heating, sunlight etc.)!

- The aluminum frames of the orthosis are coated and can be cleaned using pH-neutral soap. We recommend washing the cushions and straps by hand at 30 °C.

This product has been tested by our in-house quality management system. However, should you have any complaints, please contact your local medical supply retailer.

#### Part of the body this product is used for

In accordance with the indications (knee). Please see intended use.

#### Assembly and fitting instructions

The SecuTec Genu is supplied as a pre-assembled orthosis in standard sizes. The first fitting, which includes the hinge adjustments and the static fitting, must be carried out by a trained specialist.

#### Maintenance instructions

If handled and cared for correctly, the product is practically maintenance-free.

#### Technical specifications / parameters

The SecuTec Genu is a functional orthosis for the knee joint. It consists of a frame, hinges, straps, quick-fit fasteners and cushions.

#### Notes on reuse

This product is intended for your personal use only.

#### Warranty

The statutory regulations of the country of purchase apply. Please first contact the retailer from whom you obtained the product directly in the event of a potential claim under the warranty. The product must be cleaned before submitting it in case of warranty claims. Warranty may be limited or excluded if

the instructions on how to handle and care for the SecuTec Genu have not been observed. Warranty is excluded if:

- The product was not used according to the indication
- The instructions given by the specialist were not observed
- The product was modified arbitrarily

#### Duty to report

Due to regional legal regulations, you are obliged to immediately report any serious incident involving the use of this medical product to both the manufacturer and the relevant authorities. Our contact details can be found on the back of this brochure.

#### Disposal

Once you have finished using the product, please dispose of it according to local specifications.

#### Material content

Aluminum, Polyamide (PA), Polyurethane (PUR), Polyoxymethylene (POM), Polyethylene terephthalate (PET), Acryl Isocyanate (Powder paint), High Grade Steel, Polyamide glass fiber reinforced (PA-GF), Ethylene Vinyl Acetate (EVA), Polyester (PES), Acrylate, Brass (CuZn), Iron (Fe), Synthetic resin, Plastic film

 – Medical Device

 – DataMatrix Unique Device Identifier UDI

<sup>1</sup> Orthosis = orthopedic aid used to stabilize, relieve, immobilize, control, or correct the limbs or the torso

<sup>2</sup> A specialist is any person who is authorized, according to their state's regulations, to fit and provide instructions on the use of supports and orthoses.

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